


# PRODUCT DATA SHEET

Revision 1	Inoculating Needle	
Validfor Item-No.:	6030000/6030001	

1.	Description / Specification	
1.1	Description	Inoculating Needle 6030000: Inoculating Needle, 1µl, 221 mm in length, white, sterile 6030001: Inoculating Needle, 1µl, 221 mm in length, white, non-sterile
1.2	Dimensions	see Customer Drawing
1.3	Volume	NA
1.4	Material / Resin	PP (Polypropylene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	white
1.6	Sterilization	6030000: sterile 6030001: non-sterile
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for bacteriology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +121 °C
2.3	Autoclavability	YES
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-
3.	Packaging	
3.1	Pieces / Bag	25
3.2	Pieces / Case	1000
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
4.1	Laboratory use only	

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

