


# PRODUCT DATA SHEET

Revision 1	10 $\mu$ L Extended Pipet tips	
Validfor Item-No.:	3015000/3015010/3015001/3015011/3015100/3015110/ 3015101/3015111	

1.	Description / Specification	
1.1	Description	10 $\mu$ L Extended Pipet tips, 0.5 ~ 10 $\mu$ L, clear  3015000 10 $\mu$ L Extended Pipet tips,universal,, clear, non-sterile, bagged 3015010 10 $\mu$ L Extended Pipet tips,universal,, clear, filter, non-sterile, bagged 3015001 10 $\mu$ L Extended Pipet tips,universal,, clear, sterile, boxed 3015011 10 $\mu$ L Extended Pipet tips,universal,, clear, filter, sterile, boxed 3015100 10 $\mu$ L Extended Pipet tips, low retention, clear, non-sterile, bagged 3015110 10 $\mu$ L Extended Pipet tips, low retention, clear, filter, non-sterile, bagged 3015101 10 $\mu$ L Extended Pipet tips, low retention, clear, sterile, boxed 3015111 10 $\mu$ L Extended Pipet tips, low retention, clear, filter, sterile, boxed
1.2	Dimensions	Flask: see Customer Drawing
1.3	Volume	Work volume: 0.5 ~ 10 $\mu$ L
1.4	Material / Resin	Pipette tip: PP (Polypropylene), free of heavy metal Filter: PE (Polyethylene) Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: clear Filter: white
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +121 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	3015000/3015100	3015010/3015110	3015001/3015011/ 3015101/3015111
3.1	Pieces / Bag	1000	1000	96
3.2	Pieces / Box	20	10	50
3.3	Lot-No.	DDD YY XXX (day, year, product No.)		
3.4	Other Information	Certificate of Quality to download		

4.	Other Information	
4.1	Laboratory use only	

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.